

# Turner Lab Notebook

Callum Malcolm

## Tuesday 02-05-2023

1. Thawed N2, RBL1, and RAMOS into T75 flasks  
Note: N2 and RBL1 flasks seeded with feeders

## Wednesday 03-05-2023

1. Splitting cells

Cell Line	Cells/mL	Cells/mL Seeded	CS volume	Final Volume	Last day split	P number
N2	$8.93 \times 10^5$	$3 \times 10^5$ .	6.71 mL	12 mL	Tues 02-05-23	2
RBL1	$7.73 \times 10^5$	$3 \times 10^5$ .	7.76 mL	12 mL	Tues 02-05-23	2


2. Expanded N2-Barcodes to T175
3. Seeded feeders into 96-well
  - a. Diluted feeders to  $3.2 \times 10^4$  cells/mL
  - b. Add  $100 \mu\text{L}$  cell suspension to each well

## Monday 19-06-2023

1. Splitting cells
  - Split N2/RBL1/JIRE
  - Added 10 mL media to BLLW
    - Lots of dead cells observed
2. EC50 - ID:
  - Cylophosphamide EC50 for N2 and RBL1

## 96 Well Plate Map

EC50\_230619\_plate1



### Drug Dilutions

Drug	Dilution	Concentration	Units	Starting_Concentration
Cyclophosphamide	1	20	uM	20000
Cyclophosphamide	2	2	uM	2000
Cyclophosphamide	3	0.2	uM	200
Cyclophosphamide	4	0.02	uM	20
Cyclophosphamide	5	0.002	uM	2
Cyclophosphamide	6	0.0002	uM	0.2
Cyclophosphamide	7	2e-05	uM	0.02
Cyclophosphamide	8	2e-06	uM	0.002
Cyclophosphamide	9	2e-07	uM	0.0002

**Tuesday 20-06-2023**

### RNA Extraction for SP Test 3 Sample

- Sample collected:
- Low number of cells so used Pico-Pure Kit
- 1 sample (n = 1)

*Note before starting:*

- Pre-heat block to 42C
- Pre-make DNase treatment master mix (1.1n)

Component	Volume	MM volume
DNase I Stock	5 $\mu$ L	5.1 $\mu$ L
DNA Digestion Buffer	35 $\mu$ L	38.5 $\mu$ L

## Protocol

1. Extract cells with 100  $\mu$ L of Extraction Buffer (XB)
  - Resuspend the cell pellet gently by pipetting
  - DO NOT VORTEX
2. Incubate at 42C for 30min
3. Centrifuge sample at 3000xg for two minutes
4. Collect supernatant
  - Can stop here and freeze RNA at 70C
5. Pre-condition the RNA Purification Column
  - Add 250  $\mu$ L Conditioning Buffer onto the purification column filter membrane
  - Incubate RNA Purification Column with CB for 5min at RT
  - Centrifuge purification column at 16000xg for 1 min
6. Pipette 100  $\mu$ L of 70% Ethanol (EtOH) into cell extract
  - Mix well by pipetting
  - DO NOT VORTEX
7. Add cell extract to column
8. Centrifuge column for 2 minutes at 100 x g and then 16000 x g for 30s
9. Pipette 100  $\mu$ L Wash Buffer 1 (W1) into the purification column and centrifuge for 1min at 8000 x g
10. Add 40  $\mu$ L of DNase treatment master mix to purification tube
  - Add directly onto membrane
11. Incubate at RT for 15 min
12. Pipette 40  $\mu$ L Wash Buffer 1 (W1) onto purification column membrane and centrifuge at 8000 x g for 15 seconds
13. Pipette 100  $\mu$ L Wash Buffer 1 (W2) into purification column membrane and centrifuge for one minute at 8000 x g
14. Pipette 100  $\mu$ L Wash Buffer 1 (W2) into purification column membrane and centrifuge for two minutes at 16000 x g
  - Check if any wash buffer remains in column and centrifuge again at 16000 x g for one minute to clear all the liquid
15. Transfer purification column to new 0.5mL tube
  - Insert purification column/0.5mL assembly into 1.5mL tube for centrifugation
16. Pipette EB directly onto membrane of the purification column
  - Gently touch tip to surface of membrane to ensure maximum absorption of EB into membrane
  - Use 11  $\mu$ L to 30  $\mu$ L EB
17. Incubate purification column for one minute at room temperature

18. Centrifuge column for one minute at 1000 x g to distribute EB in column
19. Centrifuge column for one minute at 16000 x g to elute RNA
20. Start RT immediately or store at -80C until use