

PicoPure RNA Isolation Kit

Callum Malcolm

2023-06-13

Note before starting:

- Pre-heat block to 42C
- Pre-make DNase treatment master mix (1.1n)

Component	Volume	MM volume
DNase I Stock	5 μ L	
DNA Digestion Buffer	35 μ L	

Protocol

1. Extract cells with 100 μ L of Extraction Buffer (XB)
 - Resuspend the cell pellet gently by pipetting
 - DO NOT VORTEX
2. Incubate at 42C for 30min
3. Centrifuge sample at 3000xg for two minutes
4. Collect supernatant
 - Can stop here and freeze RNA at 70C
5. Pre-condition the RNA Purification Column
 - Add 250 μ L Conditioning Buffer onto the purification column filter membrane
 - Incubate RNA Purification Column with CB for 5min at RT
 - Centrifuge purification column at 16000xg for 1 min
6. Pipette 100 μ L of 70% Ethanol (EtOH) into cell extract
 - Mix well by pipetting
 - DO NOT VORTEX
7. Add cell extract to column
8. Centrifuge column for 2 minutes at 100 x g and then 16000 x g for 30s
9. Pipette 100 μ L Wash Buffer 1 (W1) into the purification column and centrifuge for 1min at 8000 x g
10. Add 40 μ L of DNase treatment master mix to purification tube
 - Add directly onto membrane
11. Incubate at RT for 15 min

12. Pipette 40 μ L Wash Buffer 1 (W1) onto purification column membrane and centrifuge at 8000 x g for 15 seconds
13. Pipette 100 μ L Wash Buffer 1 (W2) into purification column membrane and centrifuge for one minute at 8000 x g
14. Pipette 100 μ L Wash Buffer 1 (W2) into purification column membrane and centrifuge for two minutes at 16000 x g
 - Check if any wash buffer remains in column and centrifuge again at 16000 x g for one minute to clear all the liquid
15. Transfer purification column to new 0.5mL tube
 - Insert purification column/0.5mL assembly into 1.5mL tube for centrifugation
16. Pipette EB directly onto membrane of the purification column
 - Gently touch tip to surface of membrane to ensure maximum absorption of EB into membrane
 - Use 11 μ L to 30 μ L EB
17. Incubate purification column for one minute at room temperature
18. Centrifuge column for one minute at 1000 x g to distribute EB in column
19. Centrifuge column for one minute at 16000 x g to elute RNA
20. Start RT immediately or store at -80C until use