

RT-PCR WILDseq prep Protocol

Callum Malcolm

2023-06-09

RNA Extraction

Notes before starting

Use *Quick*-RNA Microprep Kit (cat# R1050/R1051)

- Use Pico-Pure for small cell pellets
- Check that all kit components have been properly constituted
- Discard flow through after every spin
- Spin columns for 1 min at max rpm (13.2rpm)
- Prepare DNase mix
 - 5 μ L DNase I + 35 μ L DNA Digestion Buffer per sample
 - Make master mix at 1.1n

Protocol:

1. Add RNA Lysis buffer
 - For cell amount 10^5 use 100 μ L
 - For cell amount 10^6 use 300 μ L
2. Add 1 volume of 100% Ethanol and mix well
 - Pipette up and down
 - Vortex for <30sec
 - Repeat until pellet is fully dissolved
3. Transfer mix into spin column
4. Wash column with 400 μ L RNA wash buffer
5. Add 40 μ L DNase I to spin column and incubate at RT for 30C for 15min
6. Add 400 μ L RNA Prep Buffer and spin
7. Add 700 μ L RNA Wash Buffer and spin
8. Add 400 μ L RNA Wash buffer and spin
9. Spin column to dry and transfer to nuclease-free tube
10. 15 μ L DNase/RNase-Free water to column matrix and centrifuge

Reverse Transcription

Notes Before Starting:

- Label PCR strip tubes
- Thaw RT reagents
- Uses WILDseq specific RT primer
- Dilute RT primer to $2\mu\text{M}$
 - Do this in a PCR hood
 - **IMPORTANT** Do not mix primers and carefully note which RT primer is added to which sample
- Make RT and ExoI master mixes (MM) before starting
 - MM amounts are 1.1n

Protocol:

Note: Use “rt” protocol in callum folder on thermocycler for heating steps

1. In PCR strip tubes dilute 1-5ug RNA in a total volume of $10\mu\text{L}$ of RNase/DNase-free water
2. Add 1 ul of RT primer
3. Add 1 ul of dNTPs
4. Denature at 65 C for 5 mins in PCR machine
5. Add straight to ice for 2 mins
6. Spin to get liquid to bottom of tube
7. Prepare master mix of RT enzyme and buffers (Make prior to starting)

Component	Volume	MM volume
5x SSIV Buffer	$4\mu\text{L}$	
SSIV RT	$1\mu\text{L}$	
100 mM DTT	$1\mu\text{L}$	
RNase Out	$1\mu\text{L}$	

8. Add $7\mu\text{L}$ of RT MM prepared above to each sample and mix
 - Spin briefly to get liquid to bottom of the tube
9. In PCR Machine: incubate at 53C for 10 mins followed by 80C for 10 mins
10. Add $3\mu\text{L}$ of ExoI MM to each sample

Component	Volume	MM volume
Thermolabile Exonuclease I	$1\mu\text{L}$	
NEBuffer r3.1*	$2\mu\text{L}$	

* Most PCR buffers are compatible

11. In PCR Machine: heat at 37C for 4 mins followed by 80C for 1 min
 - Spin briefly to get liquid to bottom of the tube
12. Add $1\mu\text{L}$ of RNase H and incubate at 37C for 20 mins
13. Dilute final cDNA 1 in 2 with DNase/RNase-free water

PCR Amplification of cDNA

1. Prepare mastermix of below reagents at 1.1x the number of samples

Component	Volume	MM volume
Q5 Master Mix	25 μ L	. \
P7_rev	2.5 μ L	
DNase/RNase H20	13.5 μ L	
Q5 polymerase	0.5 μ L	
Total Master Mix volume	40 μL	
cDNA from above	10 μ L	

2. Perform PCR using the following parameters:

Steps	Time
Step 1: 98C	30s
Step 2: 98C	10s
Step 3: 61C	30s
Step 4: 72C	30s
Step 5: 72C	2 min
Step 6: 12C	Hold

20-25 cycles of steps 2-4

3. Check each PCR reaction on TapeStation
 - Amplicon should be ~200bp