

# RT-PCR WILDseq prep Protocol

Callum Malcolm

2023-06-09

## RNA Extraction

### Notes before starting

Use *Quick*-RNA Microprep Kit (cat# R1050/R1051)

- Use Pico-Pure for small cell pellets
- Check that all kit components have been properly constituted
- Discard flow through after every spin
- Spin columns for 1 min at max rpm (13.2rpm)
- Prepare DNase mix
  - 5  $\mu$ L DNase I + 35  $\mu$ L DNA Digestion Buffer per sample
  - Make master mix at 1.1n

### Protocol:

1. Add RNA Lysis buffer
  - For cell amount  $10^5$  use 100  $\mu$ L
  - For cell amount  $10^6$  use 300  $\mu$ L
2. Add 1 volume of 100% Ethanol and mix well
  - Pipette up and down
  - Vortex for <30sec
  - Repeat until pellet is fully dissolved
3. Transfer mix into spin column
4. Wash column with 400  $\mu$ L RNA wash buffer
5. Add 40  $\mu$ L DNase I to spin column and incubate at RT for 30C for 15min
6. Add 400  $\mu$ L RNA Prep Buffer and spin
7. Add 700  $\mu$ L RNA Wash Buffer and spin
8. Add 400  $\mu$ L RNA Wash buffer and spin
9. Spin column to dry and transfer to nuclease-free tube
10. 15  $\mu$ L DNase/RNase-Free water to column matrix and centrifuge

# Reverse Transcription

## Notes Before Starting:

- Label PCR strip tubes
- Thaw RT reagents
- Uses WILDseq specific RT primer
- Dilute RT primer to 2 $\mu$ M
  - Do this in a PCR hood
  - **IMPORTANT** Do not mix primers and carefully note which RT primer is added to which sample
- Make RT and ExoI master mixes (MM) before starting
  - MM amounts are 1.1n

## Protocol:

*Note:* Use “rt” protocol in callum folder on thermocycler for heating steps

1. In PCR strip tubes dilute 1-5ug RNA in a total volume of 10  $\mu$ L of RNase/DNase-free water
2. Add 1 ul of RT primer
3. Add 1 ul of dNTPs
4. Denature at 65 C for 5 mins in PCR machine
5. Add straight to ice for 2 mins
6. Spin to get liquid to bottom of tube
7. Prepare master mix of RT enzyme and buffers (Make prior to starting)

Component	Volume	MM volume
5x SSIV Buffer	4 $\mu$ L	
SSIV RT	1 $\mu$ L	
100 mM DTT	1 $\mu$ L	
RNase Out	1 $\mu$ L	

8. Add 7  $\mu$ L of RT MM prepared above to each sample and mix
  - Spin briefly to get liquid to bottom of the tube
9. In PCR Machine: incubate at 53C for 10 mins followed by 80C for 10 mins
10. Add 3  $\mu$ L of ExoI MM to each sample

Component	Volume	MM volume
Thermolabile Exonuclease I	1 $\mu$ L	
NEBuffer r3.1*	2 $\mu$ L	

\* Most PCR buffers are compatible

11. In PCR Machine: heat at 37C for 4 mins followed by 80C for 1 min
  - Spin briefly to get liquid to bottom of the tube
12. Add 1  $\mu$ L of RNase H and incubate at 37C for 20 mins
13. Dilute final cDNA 1 in 2 with DNase/RNase-free water

## PCR Amplification of cDNA

1. Prepare mastermix of below reagents at 1.1x the number of samples

Component	Volume
5x Q5 Reaction Buffer	10 $\mu$ L
GC Enhancer	10 $\mu$ L
dNTPs (10mM)	1 $\mu$ L
P5_R1_fwd (100 $\mu$ M)	2.5 $\mu$ L
P7_rev	2.5 $\mu$ L
DNAse/RNAse H20	13.5 $\mu$ L
Q5 polymerase	0.5 $\mu$ L
<b>Total Master Mix volume</b>	<b>40 <math>\mu</math>L</b>
cDNA from above	10 $\mu$ L

2. Perform PCR using the following parameters:

Steps	Time
Step 1: 98C	30s
<b>Step 2: 98C</b>	<b>10s</b>
<b>Step 3: 61C</b>	<b>30s</b>
<b>Step 4: 72C</b>	<b>30s</b>
Step 5: 72C	2 min
Step 6: 12C	Hold

**20-25 cycles of steps 2-4**

3. Run out 5  $\mu$ L of each PCR reaction on 1.5% agarose gel

- Amplicon should be ~200bp

## Running PCR Gel protocol

- 1.