

2023-05-22 N2-Barcode Drug Pressure

Drug Diltutions

Make 1 μ M working stock of each drug

* Doxorubicin: 1 μ L 10mM stock in 10mL media

* Methotrexate: 1 μ L 10mM stock in 10mL media

* Vincristine: 333.33 μ L 30mM stock in 10mL media

Doxorubicin

Drug	Drug ID	Dilution volume	Media Volume	[Final Drug]
Doxorubicin	Dox 1	3600 μ L (from WS)	8400 μ L	300nm
Doxorubicin	Dox 2	6000 μ L (from Dox 1)	6000 μ L	150nm
Doxorubicin	Dox 3	6000 μ L (from Dox 2)	6000 μ L	75nm
Doxorubicin	Dox 4	6000 μ L (from Dox 3)	6000 μ L	37.5nm
Doxorubicin	Dox 5	6000 μ L (from Dox 4)	6000 μ L	18.75nm

Methotrexate

Drug	Drug ID	Dilution volume	Media Volume	[Final Drug]
Methotrexate	Meth 1	1200 μ L (from WS)	10800 μ L	100nm
Methotrexate	Meth 2	6000 μ L (from Meth 1)	6000 μ L	50nm
Methotrexate	Meth 3	6000 μ L (from Meth 2)	6000 μ L	20nm
Methotrexate	Meth 4	6000 μ L (from Meth 3)	6000 μ L	10nm
Methotrexate	Meth 5	6000 μ L (from Meth 4)	6000 μ L	5nm

Vincristine

Drug	Drug ID	Dilution Volume	Media Volume	[Final Drug]
Vincristine	Vin 1	48 μ L (from WS)	11952 μ L	4nm
Vincristine	Vin 2	6000 μ L (from Vin 1)	6000 μ L	2nm
Vincristine	Vin 3	6000 μ L (from Vin 2)	6000 μ L	1nm
Vincristine	Vin 4	6000 μ L (from Vin 3)	6000 μ L	0.5nm
Vincristine	Vin 5	6000 μ L (from Vin 4)	6000 μ L	0.25nm

Instructions

1. Make drug dilutions and DMSO tubes (3 μ L DMSO in 30mL)
2. Collect wells into labeled 15mL tubes
3. Count cells
4. Spin down volume of CS which $\leq 1.0 \times 10^6$
5. Resuspend in associated drug mix and add back to well