

# Turner Lab Notebook

Callum Malcolm

## Contents

<b>Tuesday 02-05-2023</b>	<b>1</b>
<b>Wednesday 03-05-2023</b>	<b>1</b>
<b>Monday 19-06-2023</b>	<b>2</b>
Seed EC50 . . . . .	2
<b>Tuesday 20-06-2023</b>	<b>3</b>
RNA Extraction for SP Test 3 Sample . . . . .	3
Tapestation_230620 . . . . .	4
Reverse Transcription . . . . .	4
PCR Amplification of cDNA . . . . .	5
<b>Wednesday 21-06-2023</b>	<b>6</b>
Tapestation_230619_SPTTest3_cDNA . . . . .	6
<b>Thursday 22-06-2023</b>	<b>6</b>
Collect EC50 Plate . . . . .	6
<b>Friday 23-06-23</b>	<b>6</b>
Split Cells . . . . .	6
Seed EC50 . . . . .	7

## Tuesday 02-05-2023

1. Thawed N2, RBL1, and RAMOS into T75 flasks  
Note: N2 and RBL1 flasks seeded with feeders

## Wednesday 03-05-2023

1. Splitting cells

Cell Line	Cells/mL	Cells/mL Seeded	CS volume	Final Volume	Last day split	P number
N2	8.93x10 <sup>5</sup>	3x10 <sup>5</sup> .	6.71 mL	12 mL	Tues 02-05-23	2
RBL1	7.73x10 <sup>5</sup>	3x10 <sup>5</sup> .	7.76 mL	12 mL	Tues 02-05-23	2

2. Expanded N2-Barcodes to T175
3. Seeded feeders into 96-well
  - a. Diluted feeders to 3.2x10<sup>4</sup> cells/mL
  - b. Add 100μL cell suspension to each well

Monday 19-06-2023

1. Splitting cells


- Split N2/RBL1/JIRE
- Added 10 mL media to BLLW
  - Lots of dead cells observed

Seed EC50

- Plate ID: EC50\_230619\_Plate1
- Cyclophosphamide (Stock: 20mM) EC50 for N2 and RBL1

96 Well Plate Map

EC50\_230619\_plate1



Drug Dilutions

Drug	Dilution	Concentration	Units	Starting Concentration
Cyclophosphamide	1	20	uM	20000
Cyclophosphamide	2	2	uM	2000
Cyclophosphamide	3	0.2	uM	200
Cyclophosphamide	4	0.02	uM	20
Cyclophosphamide	5	0.002	uM	2
Cyclophosphamide	6	0.0002	uM	0.2
Cyclophosphamide	7	2e-05	uM	0.02
Cyclophosphamide	8	2e-06	uM	0.002
Cyclophosphamide	9	2e-07	uM	0.0002

**Tuesday 20-06-2023**

### **RNA Extraction for SP Test 3 Sample**

- Sample collected:
- Low number of cells so used Pico-Pure Kit
- 1 sample (n = 1)

*Note before starting:*

- Pre-heat block to 42C
- Pre-make DNase treatment master mix (1.1n)

Component	Volume	MM volume
DNase I Stock	5 $\mu$ L	5.5 $\mu$ L
DNA Disgestion Buffer	35 $\mu$ L	38.5 $\mu$ L

### **RNA Extraction Pico-Pure Protocol**

1. Extract cells with 100  $\mu$ L of Extraction Buffer (XB)
  - Resuspend the cell pellet gently by pipetting
  - DO NOT VORTEX
2. Incubate at 42C for 30min
3. Centrifuge sample at 3000xg for two minutes
4. Collect supernatant
  - Can stop here and freeze RNA at 70C
5. Pre-condition the RNA Purification Column
  - Add 250  $\mu$ L Conditioning Buffer onto the purification column filter membrane
  - Incubate RNA Purification Column with CB for 5min at RT
  - Centrifuge purification column at 16000xg for 1 min
6. Pipette 100  $\mu$ L of 70% Ethanol (EtOH) into cell extract
  - Mix well by pipetting
  - DO NOT VORTEX
7. Add cell extract to column
8. Centrifuge column for 2 minutes at 100 x g and then 16000 x g for 30s
9. Pipette 100  $\mu$ L Wash Buffer 1 (W1) into the purification column and centrifuge for 1min at 8000 x g
10. Add 40  $\mu$ L of DNase treatment master mix to purification tube
  - Add directly onto membrane
11. Incubate at RT for 15 min
12. Pipette 40  $\mu$ L Wash Buffer 1 (W1) onto purification column membrane and centrifuge at 8000 x g for 15 seconds

13. Pipette 100  $\mu$ L Wash Buffer 1 (W2) into purification column membrane and centrifuge for one minute at 8000 x g
14. Pipette 100  $\mu$ L Wash Buffer 1 (W2) into purification column membrane and centrifuge for two minutes at 16000 x g
  - Check if any wash buffer remains in column and centrifuge again at 16000 x g for one minute to clear all the liquid
15. Transfer purification column to new 0.5mL tube
  - Insert purification column/0.5mL assembly into 1.5mL tube for centrifugation
16. Pipette EB directly onto membrane of the purification column
  - Gently touch tip to surface of membrane to ensure maximum absorption of EB into membrane
  - Use 11  $\mu$ L to 30  $\mu$ L EB
17. Incubate purification column for one minute at room temperature
18. Centrifuge column for one minute at 1000 x g to distribute EB in column
19. Centrifuge column for one minute at 16000 x g to elute RNA
20. Started PCR and completed Tapestation

## Tapestation\_\_230620

### High-Sensitivity RNA ScreenTape Protocol

#### Results

- Inconclusive, proceeded with PCR

## Reverse Transcription

### Notes Before Starting:

- Label PCR strip tubes
- Thaw RT reagents
- Uses WILDseq specific RT primer
- Dilute RT primer to 2 $\mu$ M
  - Do this in a PCR hood
  - **IMPORTANT** Do not mix primers and carefully note which RT primer is added to which sample
- Make RT and ExoI master mixes (MM) before starting
  - MM amounts are 1.1n

### Protocol:

*Note:* Used “rt” protocol in callum folder on thermocycler for heating steps

1. In PCR strip tubes diluted 1-5ug RNA in a total volume of 10  $\mu$ L of RNase/DNase-free water
2. Added 1 ul of RT primer
3. Added 1 ul of dNTPs
4. Denatured at 65 C for 5 mins in PCR machine
5. Added straight to ice for 2 mins

6. Spun to get liquid to bottom of tube
7. Prepared master mix of RT enzyme and buffers (Made prior to starting)

Component	Volume	MM volume
5x SSIV Buffer	4 $\mu$ L	4.4 $\mu$ L
SSIV RT	1 $\mu$ L	1.1 $\mu$ L
100 mM DTT	1 $\mu$ L	1.1 $\mu$ L
RNase Out	1 $\mu$ L	1.1 $\mu$ L

8. Added 7  $\mu$ L of RT MM prepared above to each sample and mix
  - Spun briefly to get liquid to bottom of the tube
9. In PCR Machine: incubated at 53C for 10 mins followed by 80C for 10 mins
10. Added 3  $\mu$ L of ExoI MM to each sample

Component	Volume	MM volume
Thermolabile Exonuclease I	1 $\mu$ L	1.1 $\mu$ L
NEBuffer r3.1*	2 $\mu$ L	2.2 $\mu$ L

- Most PCR buffers are compatible
11. In PCR Machine: heat at 37C for 4 mins followed by 80C for 1 min
    - Spin briefly to get liquid to bottom of the tube
  12. Add 1  $\mu$ L of RNase H and incubate at 37C for 20 mins
  13. Dilute final cDNA 1 in 2 with DNase/RNase-free water

## PCR Amplification of cDNA

1. Prepare mastermix of below reagents at 1.1x the number of samples

Component	Volume	MM volume
Q5 Master Mix	25 $\mu$ L	
P5_R1_fwd (100 $\mu$ M)	2.5 $\mu$ L	
P7_rev	2.5 $\mu$ L	
DNase/RNase H20	13.5 $\mu$ L	
Q5 polymerase	0.5 $\mu$ L	
<b>Total Master Mix volume</b>	<b>40 <math>\mu</math>L</b>	
cDNA from above	10 $\mu$ L	

2. Perform PCR using the following parameters:

Steps	Time
Step 1: 98C	30s
<b>Step 2: 98C</b>	<b>10s</b>
<b>Step 3: 61C</b>	<b>30s</b>
<b>Step 4: 72C</b>	<b>30s</b>

Steps	Time
Step 5: 72C	2 min
Step 6: 12C	Hold

**20-25 cycles of steps 2-4**

Ran overnight and left at 12C

**Wednesday 21-06-2023**

**Tapestation\_230619\_SPTTest3\_cDNA**

[High-Sensitivity DNA ScreenTape Protocol](#)

[Results](#)

\* Did not work, bands present but faint and not at 200bp


\* Potential solutions:

- Repeat PCR with increased cycles
- Repeat RNA extration with trizol protocol (Anna)

**Thursday 22-06-2023**

**Collect EC50 Plate**

- Plate ID: EC50\_230619\_Plate1




2. Split JIRE 1/2
3. Split N2/RBL1 1/4

#### Seed EC50

- EC50 Plate ID: EC50\_230623\_Plate1

#### 96 Well Plate Map

EC50\_230623\_plate1



#### Drug Dilutions

Drug	Dilution	Concentration	Units	Starting_Concentration
Cyclophosphamide	1	20	uM	20000
Cyclophosphamide	2	13.3	uM	13300
Cyclophosphamide	3	8.89	uM	8890
Cyclophosphamide	4	5.93	uM	5930
Cyclophosphamide	5	3.95	uM	3950
Cyclophosphamide	6	2.63	uM	2630
Cyclophosphamide	7	1.76	uM	1760
Cyclophosphamide	8	1.17	uM	1170
Cyclophosphamide	9	0.78	uM	780