

# Turner Lab Notebook

Callum Malcolm

## Contents

<b>MARCH 2023</b>	<b>2</b>
<b>Wednesday 15-03-2023</b>	<b>2</b>
Cell Culture . . . . .	2
N2 Growth Curve . . . . .	2
<b>Thursday 16-03-2023</b>	<b>2</b>
N2/BLLW Growth Curve Day 1 . . . . .	2
<b>Friday 17-03-2023</b>	<b>3</b>
N2/BLLW Growth Curve Day 2 . . . . .	3
Split Cells . . . . .	3
<b>Saturday 18-03-2023</b>	<b>3</b>
N2/BLLW Growth Curve Day 3 . . . . .	3
<b>Sunday 18-03-2023</b>	<b>4</b>
N2/BLLW Growth Curve Day 4 . . . . .	4
<b>Monday 20-03-2023</b>	<b>4</b>
N2 Growth Curve Day 5 . . . . .	4
Split Cells . . . . .	4
Seed BLLW . . . . .	4
<b>Tuesday 21-03-2023</b>	<b>4</b>
<b>Friday 24-03-2023</b>	<b>4</b>
Seed Cells - N2/RBL1/JIRE . . . . .	4
<b>Monday 27-03-2023</b>	<b>4</b>
Split Cells . . . . .	4
<b>Wednesday 29-03-2023</b>	<b>5</b>
Split Cells . . . . .	5
<b>Friday 31-03-2023</b>	<b>5</b>
Split Cells . . . . .	5
<b>APRIL 2023</b>	<b>5</b>
<b>Monday 03-04-2023</b>	<b>5</b>
Seed Cells - N2/RBL1/JIRE . . . . .	5
Seed EC50 . . . . .	5

<b>Wednesday 05-04-2023</b>	<b>5</b>
Seed Cells - N2/RBL1/JIRE . . . . .	5
<b>Thursday 06-04-2023</b>	<b>5</b>
Collect EC50 Plate . . . . .	5
<b>Tuesday 11-04-2023</b>	<b>6</b>
Seed EC50_230411 . . . . .	6
<b>Wednesday 12-04-2023</b>	<b>7</b>
Cell Culture . . . . .	7
Side Population - N2_230412 . . . . .	7
<b>Tuesday 02-05-2023</b>	<b>7</b>
Seeding Cells - N2/RBL1/RAMOS . . . . .	7
<b>Wednesday 03-05-2023</b>	<b>7</b>
Split Cells . . . . .	7
<b>Thursday 04-05-2023</b>	<b>7</b>
<b>Monday 19-06-2023</b>	<b>7</b>
Split Cells . . . . .	7
Seed EC50 . . . . .	8
<b>Saturday 20-05-2023</b>	<b>8</b>
<b>JUNE 2023</b>	<b>8</b>
<b>Tuesday 20-06-2023</b>	<b>8</b>
RNA Extraction for SP Test 3 Sample . . . . .	8
Tapestation_230620 . . . . .	10
Reverse Transcription . . . . .	10
PCR Amplification of cDNA . . . . .	11
<b>Wednesday 21-06-2023</b>	<b>12</b>
Tapestation_230619_SPTTest3_cDNA . . . . .	12
<b>Thursday 22-06-2023</b>	<b>12</b>
Collect EC50 Plate . . . . .	12
<b>Friday 23-06-23</b>	<b>12</b>
Split Cells . . . . .	12
Seed EC50 . . . . .	12
<b>Monday 26-06-23</b>	<b>13</b>
Collect EC50 Plate . . . . .	13
<b>JULY 2023</b>	<b>14</b>
<b>Monday 24-07-2023</b>	<b>14</b>
Seeding Cells - RAMOS/FDC . . . . .	14
Made New Media: AR-1 . . . . .	14
<b>Tuesday 25-07-2023</b>	<b>15</b>
Seed Cells - N2-BC4/RBL1/BLLW/RAMOS-BC1 . . . . .	15

<b>Wednesday 26-07-2023</b>	<b>15</b>
Cell Culture . . . . .	15
<b>Thursday 27-07-2023</b>	<b>15</b>
Split Cells - RAMOS . . . . .	15
Lucy Drug Randomizer . . . . .	16
<b>Friday 28-07-2023</b>	<b>16</b>
<b>OCTOBER 2023</b>	<b>16</b>
<b>Friday 06-10-2023</b>	<b>16</b>
Cell Culture . . . . .	16
Made New Media: AR-2 . . . . .	16
<b>Monday 09-10-2023</b>	<b>16</b>
Cell Culture . . . . .	16
Seeded A20 . . . . .	16
Seeded FDC . . . . .	17
<b>Tuesday 10-10-2023</b>	<b>17</b>
Cell Culture . . . . .	17
Seeded N4 . . . . .	17
Seeded RAMOS-BC . . . . .	17
Checked A20 . . . . .	17
<b>Wednesday 11-10-2023</b>	<b>17</b>
Cell Culture . . . . .	17
A20 Culture . . . . .	17
RAMOS-BC4 Culture . . . . .	17
N4 Culture . . . . .	17
<b>Thursday 12-10-2023</b>	<b>18</b>
Cell Culture . . . . .	18
Ramos-BC4 . . . . .	18
A20 . . . . .	18
N4 . . . . .	18
<b>Friday 13-10-2023</b>	<b>18</b>
Cell Culture . . . . .	18
A20 . . . . .	18
EC50_231013 . . . . .	18
EC50 Plate ID: EC50_231013_Plate1 . . . . .	18
EC50 Plate ID: EC50_231013_Plate2 . . . . .	19
<b>Monday 16-10-2023</b>	<b>20</b>
Cell Culture . . . . .	20
A20 . . . . .	20
N4 . . . . .	20
Collect EC50_231013 . . . . .	20
Plate ID: EC50_231013_Plate1 . . . . .	20
Plate ID: EC50_231013_Plate2 . . . . .	22
<b>Thursday 19-10-2023</b>	<b>22</b>
Cell Culture . . . . .	22
A20 . . . . .	22

N4 . . . . .	22
<b>Friday 20-10-2023</b>	<b>22</b>
Cell Culture . . . . .	22
A20 . . . . .	22
N4 . . . . .	22
<b>Sunday 22-10-2023</b>	<b>23</b>
Cell Culture . . . . .	23
RAMOS Barcode . . . . .	23
<b>Monday 23-10-2023</b>	<b>23</b>
Cell Culture . . . . .	23
N4 . . . . .	23
A20 . . . . .	23
Ramos-BC 1 . . . . .	23
Ramos-BC 5 . . . . .	23
Ramos-BC 6 . . . . .	23
Seed Feeders (FDC) . . . . .	23
<b>Tuesday 23-10-2023</b>	<b>23</b>
Cell Culture . . . . .	23
N4 . . . . .	23
A20 . . . . .	24
Ramos-BC 1 . . . . .	24
Ramos-BC 5 . . . . .	24
Ramos-BC 6 . . . . .	24
FDC . . . . .	24
Transduce A20/N4 with WILDseq Library . . . . .	24
<b>Wednesday 25-10-23</b>	<b>24</b>
Attempt 1: A20/N4 Barcoding - Media change . . . . .	24
Made New Media: AR-3 . . . . .	25
<b>Thursday 26-10-23</b>	<b>25</b>
Attempt 1: A20/N4 Barcoding - Sort . . . . .	25
<b>Friday 27-10-23</b>	<b>25</b>
Cell Culture . . . . .	25
A20 . . . . .	25
N4 . . . . .	25
Ramos-BC 1 . . . . .	25
Ramos-BC 5 . . . . .	25
Ramos-BC 6 . . . . .	25
EC50_231027 . . . . .	26
Overview . . . . .	26
EC50 Protocol: . . . . .	26
EC50 Plate ID: EC50_231027_Plate1 . . . . .	26
EC50 Plate ID: EC50_231027_Plate2 . . . . .	27
<b>Monday 27-10-23</b>	<b>28</b>
Cell Culture . . . . .	28
A20 . . . . .	28
N4 . . . . .	28
Ramos-BC 1 . . . . .	28

Ramos-BC 5 . . . . .	29
Ramos-BC 6 . . . . .	29
Seed Feeders . . . . .	29
Collect EC50_231027 . . . . .	29
Plate ID: EC50_231027_Plate . . . . .	29
<b>Tuesday 31-10-2023</b>	<b>29</b>
Attempt 2: Transduce A20/N4 with WILDseq Library . . . . .	29
Drug Pressure Experiment 1: RAMOS Day 1 . . . . .	30
Overview . . . . .	30
Day 1 . . . . .	30
Seeded RAMOS BC . . . . .	30
Collected RAMOS Baseline . . . . .	31
<b>NOVEMBER 2023</b>	<b>31</b>
<b>Wednesday 01-11-2023</b>	<b>31</b>
Drug Pressure Experiment 1: RAMOS Day 2 . . . . .	31
Attempt 2: Barcoding A20/N4 . . . . .	32
<b>Thursday 02-11-2023</b>	<b>32</b>
Cell Culture . . . . .	32
A20 . . . . .	32
N4 . . . . .	32
Ramos-BC 1 . . . . .	32
Ramos-BC 5 . . . . .	32
Ramos-BC 6 . . . . .	32
Seeded Feeders . . . . .	32
Drug Pressure Experiment 1: RAMOS Day 3 . . . . .	32
<b>Friday 03-11-2023</b>	<b>33</b>
Cell Culture . . . . .	33
A20 . . . . .	33
N4 . . . . .	33
Ramos-BC 1 . . . . .	33
Ramos-BC 5 . . . . .	33
Ramos-BC 6 . . . . .	33
Drug Pressure Experiment 1: RAMOS Day 4 . . . . .	33
<b>Monday 06-11-2023</b>	<b>33</b>
Cell Culture . . . . .	33
A20 . . . . .	33
N4.1 . . . . .	33
N4.2 . . . . .	34
Ramos-BC 1 . . . . .	34
Ramos-BC 5 . . . . .	34
Ramos-BC 6 . . . . .	34
Seed Feeders . . . . .	34
Drug Pressure Experiment 1: RAMOS Day 7 . . . . .	34
<b>Tuesday 07/11/2023</b>	<b>34</b>
Made New Media: AR-4 . . . . .	34
Cell Culture . . . . .	34
A20 . . . . .	34
N4.1 . . . . .	35

N4.2 . . . . .	35
Drug Pressure Experiment 1: RAMOS Day 5 . . . . .	35
Attempt 3: Transduce N4 with WILDseq Library . . . . .	35
<b>Wednesday 08/11/2023</b>	<b>35</b>
Attempt 3: Barcoding N4 . . . . .	35
<b>Thursday 09/11/2023</b>	<b>35</b>
Attempt 3: Barcoding N4 . . . . .	35
A20 . . . . .	36
<b>Friday 10/11/2023</b>	<b>36</b>
Attempt 3: Barcoding N4 . . . . .	36
Cell Culture . . . . .	36
RAMOS-BC1 . . . . .	36
<b>Monday 13/11/2023</b>	<b>36</b>
Start RAMOS-Drug Pressure v2 . . . . .	36
Cell Culture . . . . .	38
N4-BC Attempt 3 . . . . .	38
<b>Thursday 13/11/2023</b>	<b>38</b>
RAMOS-DPv2 . . . . .	38
<b>Friday 14/11/2023</b>	<b>38</b>
Cell Cutlure . . . . .	38

**MARCH 2023**

**Wednesday 15-03-2023**

### Cell Culture


- Given N2/RBL1/BLLW PDX line from Chris Steele

### N2 Growth Curve

Looking to compare N2 growth on and off feeders

1. Seeded 2mL of N2/BLLW cell suspension at  $0.25 \times 10^6$  cells/mL

## 96 Well Plate Map



**Thursday 16-03-2023**

### N2/BLLW Growth Curve Day 1

- Cells gently mixed
- 30  $\mu$ L cell suspension taken and counted
- Only counted 1 square and multiplied by 4

Cell Line	FDC Status	Count	Cell Number
N2	+ FDC	25	2000000
N2	+ FDC	26	2080000
N2	+ FDC	12	960000
N2	- FDC	15	1200000
N2	- FDC	14	1120000
N2	- FDC	6	480000
BLLW	+ FDC	6	480000
BLLW	+ FDC	6	480000
BLLW	+ FDC	5	400000
BLLW	- FDC	3	240000
BLLW	- FDC	3	240000
BLLW	- FDC	6	480000

## Friday 17-03-2023

### N2/BLLW Growth Curve Day 2

- Cells gently mixed
- 30  $\mu$ L cell suspension taken and counted
- Only counted 1 square and multiplied by 4

Cell Line	FDC Status	Count	Cell Number
N2	+ FDC	36	2880000
N2	+ FDC	28	2240000
N2	+ FDC	38	3040000
N2	- FDC	61	4880000
N2	- FDC	38	3040000
N2	- FDC	16	1280000
BLLW	+ FDC	3	240000
BLLW	+ FDC	5	400000
BLLW	+ FDC	5	400000
BLLW	- FDC	1	80000
BLLW	- FDC	0	0
BLLW	- FDC	1	80000

### Split Cells

Cell Line	Cells/mL	Cells/mL Seeded	CS volume	Final Volume	Seeded Cells/mL
N2	$8.8 \times 10^5$	$3 \times 10^5$	6.8 mL	12 mL	$3 \times 10^5$
RBL1	$8.6 \times 10^5$	$3 \times 10^5$	11.6 mL	12 mL	$5 \times 10^5$

## Saturday 18-03-2023

### N2/BLLW Growth Curve Day 3

- Counted by Chris



## Sunday 18-03-2023

### N2/BLLW Growth Curve Day 4

- Counted by Chris

## Monday 20-03-2023

### N2 Growth Curve Day 5

- BLLW were all dead, experiment stopped
- Cells gently mixed
- 30  $\mu$ L cell suspension taken and counted
- Only counted 1 square and multiplied by 4

Cell Line	FDC Status	Count	Cell Number
N2	+ FDC	146	11680000
N2	+ FDC	112	8960000
N2	+ FDC	132	10560000
N2	- FDC	112	8960000
N2	- FDC	155	12400000
N2	- FDC	141	11280000

### Split Cells

- Split cell lines N2/RBL1/JIRE
- Dirty split 1:2

### Seed BLLW

- Given BLLW aliquot from Chris Steel
- Thawed and resuspended in 12 mL media
- Added to T75 pre-seeded with FDC cells (from Chris Steel)

## Tuesday 21-03-2023

### Large Infection in TC

- N2 Growth Curve stopped

## Friday 24-03-2023

### Seed Cells - N2/RBL1/JIRE

- Vials given by Chris Steel
- Thawed and resuspended in 12 mL media
- Added to T75 pre-seeded with FDC cells (from Chris Steel)

## Monday 27-03-2023

### Split Cells

- Split cell lines N2/RBL1/JIRE
- Dirty split 1:3

## **Wednesday 29-03-2023**

### **Split Cells**

- Split by Chris Steel.

## **Friday 31-03-2023**

### **Split Cells**

- Split cell lines N2/RBL1/JIRE
- Dirty split 1:4

## **APRIL 2023**

## **Monday 03-04-2023**

### **Large Infection in TC**

### **Seed Cells - N2/RBL1/JIRE**

- Vials given by Chris Steel
- Thawed and resuspended in 12 mL media
- Added to T75 pre-seeded with FDC cells (from Chris Steel)

### **Seed EC50**

DNW - N2/JIRE/RBL1  
- Vin/Dactolisib

## **Wednesday 05-04-2023**

### **Seed Cells - N2/RBL1/JIRE**

- Vials given by Chris Steel
- Thawed and resuspended in 12 mL media
- Added to T75 pre-seeded with FDC cells (from Chris Steel)

## **Thursday 06-04-2023**

### **Collect EC50 Plate**


- Plate ID: EC50\_230619\_Plate1  
- DNW

Tuesday 11-04-2023

Seed EC50\_230411

## 96 Well Plate Map

### EC50\_230411\_drug\_plating



#### Drug Dilutions

Drug	Dilution	Concentration	Units	Starting_Concentration
Dactolisib	1	20	uM	20000
Dactolisib	2	2	uM	2000
Dactolisib	3	0.2	uM	200
Dactolisib	4	0.02	uM	20
Dactolisib	5	0.002	uM	2
Dactolisib	6	0.0002	uM	0.2
Dactolisib	7	2e-05	uM	0.02
Dactolisib	8	2e-06	uM	0.002
Dactolisib	9	2e-07	uM	0.0002

Drug	Dilution	Concentration	Units	Starting_Concentration
Vincristine	1	30	nM	30000
Vincristine	2	3	nM	3000
Vincristine	3	0.3	nM	300
Vincristine	4	0.03	nM	30
Vincristine	5	0.003	nM	3
Vincristine	6	0.0003	nM	0.3
Vincristine	7	3e-05	nM	0.03
Vincristine	8	3e-06	nM	0.003
Vincristine	9	3e-07	nM	0.0003

- 3 plates seeded
  - Plate 1: N2

- Plate 2: JIRE
- Plate 3: RBL1

## Wednesday 12-04-2023

### Cell Culture

- Thaw RBL1

### Side Population - N2\_230412

- Hoescht
  - 2.5 / 2.75 / 3
- DNW

# MAY 2023

## Tuesday 02-05-2023

### Seeding Cells - N2/RBL1/RAMOS

1. Thawed N2, RBL1, and RAMOS into T75 flasks  
Note: N2 and RBL1 flasks seeded with feeders

## Wednesday 03-05-2023

### Split Cells

Cell Line	Cells/mL	Cells/mL Seeded	CS volume	Final Volume	Last day split	P number
N2	8.93x10 <sup>5</sup>	3x10 <sup>5</sup> .	6.71 mL	12 mL	Tues 02-05-23	2
RBL1	7.73x10 <sup>5</sup>	3x10 <sup>5</sup> .	7.76 mL	12 mL	Tues 02-05-23	2

2. Expanded N2-Barcodes to T175
3. Seeded feeders into 96-well
  - a. Diluted feeders to 3.2x10<sup>4</sup> cells/mL
  - b. Add 100μL cell suspension to each well

## Thursday 04-05-2023

## Monday 19-06-2023

### Split Cells


- Split N2/RBL1/JIRE
- Added 10 mL media to BLLW
  - Lots of dead cells observed

## Seed EC50

- Plate ID: EC50\_230619\_Plate1
- Cyclophosphamide (Stock: 20mM) EC50 for N2 and RBL1

## 96 Well Plate Map

### EC50\_230619\_plate1



## Drug Dilutions

Drug	Dilution	Concentration	Units	Starting_Concentration
Cyclophosphamide	1	20	uM	20000
Cyclophosphamide	2	2	uM	2000
Cyclophosphamide	3	0.2	uM	200
Cyclophosphamide	4	0.02	uM	20
Cyclophosphamide	5	0.002	uM	2
Cyclophosphamide	6	0.0002	uM	0.2
Cyclophosphamide	7	2e-05	uM	0.02
Cyclophosphamide	8	2e-06	uM	0.002
Cyclophosphamide	9	2e-07	uM	0.0002

Saturday 20-05-2023

JUNE 2023

Tuesday 20-06-2023

## RNA Extraction for SP Test 3 Sample

- Sample collected:

- Low number of cells so used Pico-Pure Kit
- 1 sample (n = 1)

*Note before starting:*

- Pre-heat block to 42C
- Pre-make DNase treatment master mix (1.1n)

Component	Volume	MM volume
DNase I Stock	5 $\mu$ L	5.5 $\mu$ L
DNA Digestion Buffer	35 $\mu$ L	38.5 $\mu$ L

### RNA Extraction Pico-Pure Protocol

1. Extract cells with 100  $\mu$ L of Extraction Buffer (XB)
  - Resuspend the cell pellet gently by pipetting
  - DO NOT VORTEX
2. Incubate at 42C for 30min
3. Centrifuge sample at 3000xg for two minutes
4. Collect supernatant
  - Can stop here and freeze RNA at 70C
5. Pre-condition the RNA Purification Column
  - Add 250  $\mu$ L Conditioning Buffer onto the purification column filter membrane
  - Incubate RNA Purification Column with CB for 5min at RT
  - Centrifuge purification column at 16000xg for 1 min
6. Pipette 100  $\mu$ L of 70% Ethanol (EtOH) into cell extract
  - Mix well by pipetting
  - DO NOT VORTEX
7. Add cell extract to column
8. Centrifuge column for 2 minutes at 100 x g and then 16000 x g for 30s
9. Pipette 100  $\mu$ L Wash Buffer 1 (W1) into the purification column and centrifuge for 1min at 8000 x g
10. Add 40  $\mu$ L of DNase treatment master mix to purification tube
  - Add directly onto membrane
11. Incubate at RT for 15 min
12. Pipette 40  $\mu$ L Wash Buffer 1 (W1) onto purification column membrane and centrifuge at 8000 x g for 15 seconds
13. Pipette 100  $\mu$ L Wash Buffer 1 (W2) into purification column membrane and centrifuge for one minute at 8000 x g
14. Pipette 100  $\mu$ L Wash Buffer 1 (W2) into purification column membrane and centrifuge for two minutes at 16000 x g
  - Check if any wash buffer remains in column and centrifuge again at 16000 x g for one minute to clear

all the liquid

15. Transfer purification column to new 0.5mL tube
  - Insert purification column/0.5mL assembly into 1.5mL tube for centrifugation
16. Pipette EB directly onto membrane of the purification column
  - Gently touch tip to surface of membrane to ensure maximum absorption of EB into membrane
  - Use 11  $\mu$ L to 30  $\mu$ L EB
17. Incubate purification column for one minute at room temperature
18. Centrifuge column for one minute at 1000 x g to distribute EB in column
19. Centrifuge column for one minute at 16000 x g to elute RNA
20. Started PCR and completed Tapestation

## Tapestation\_230620

### High-Sensitivity RNA ScreenTape Protocol

#### Results

- Inconclusive, proceeded with PCR

## Reverse Transcription

### Notes Before Starting:

- Label PCR strip tubes
- Thaw RT reagents
- Uses WILDseq specific RT primer
- Dilute RT primer to 2 $\mu$ M
  - Do this in a PCR hood
  - **IMPORTANT** Do not mix primers and carefully note which RT primer is added to which sample
- Make RT and ExoI master mixes (MM) before starting
  - MM amounts are 1.1n

### Protocol:

*Note:* Used “rt” protocol in callum folder on thermocycler for heating steps

1. In PCR strip tubes diluted 1-5ug RNA in a total volume of 10  $\mu$ L of RNase/DNase-free water
2. Added 1 ul of RT primer
3. Added 1 ul of dNTPs
4. Denatured at 65 C for 5 mins in PCR machine
5. Added straight to ice for 2 mins
6. Spun to get liquid to bottom of tube
7. Prepared master mix of RT enzyme and buffers (Made prior to starting)

Component	Volume	MM volume
5x SSIV Buffer	4 $\mu$ L	4.4 $\mu$ L

Component	Volume	MM volume
SSIV RT	1 $\mu$ L	1.1 $\mu$ L
100 mM DTT	1 $\mu$ L	1.1 $\mu$ L
RNAse Out	1 $\mu$ L	1.1 $\mu$ L

8. Added 7  $\mu$ L of RT MM prepared above to each sample and mix
  - Spun briefly to get liquid to bottom of the tube
9. In PCR Machine: incubated at 53C for 10 mins followed by 80C for 10 mins
10. Added 3  $\mu$ L of ExoI MM to each sample

Component	Volume	MM volume
Thermolabile Exonuclease I	1 $\mu$ L	1.1 $\mu$ L
NEBuffer r3.1*	2 $\mu$ L	2.2 $\mu$ L

- Most PCR buffers are compatible
11. In PCR Machine: heat at 37C for 4 mins followed by 80C for 1 min
    - Spin briefly to get liquid to bottom of the tube
  12. Add 1  $\mu$ L of RNAse H and incubate at 37C for 20 mins
  13. Dilute final cDNA 1 in 2 with DNase/RNAse-free water

## PCR Amplification of cDNA

1. Prepare mastermix of below reagents at 1.1x the number of samples

Component	Volume	MM volume
Q5 Master Mix	25 $\mu$ L	
P5_R1_fwd (100 $\mu$ M)	2.5 $\mu$ L	
P7_rev	2.5 $\mu$ L	
DNase/RNAse H20	13.5 $\mu$ L	
Q5 polymerase	0.5 $\mu$ L	
<b>Total Master Mix volume</b>	<b>40 <math>\mu</math>L</b>	
cDNA from above	10 $\mu$ L	

2. Perform PCR using the following parameters:

Steps	Time
Step 1: 98C	30s
<b>Step 2: 98C</b>	<b>10s</b>
<b>Step 3: 61C</b>	<b>30s</b>
<b>Step 4: 72C</b>	<b>30s</b>
Step 5: 72C	2 min
Step 6: 12C	Hold

**20-25 cycles of steps 2-4**

Ran overnight and left at 12C



**Wednesday 21-06-2023**

**Tapestation\_230619\_SPTest3\_cDNA**

High-Sensitivity DNA ScreenTape Protocol



**Results**

- \* Did not work, bands present but faint and not at 200bp
- \* Potential solutions:
  - Repeat PCR with increased cycles
  - Repeat RNA extraction with trizol protocol (Anna)

**Thursday 22-06-2023**

**Collect EC50 Plate**

- Plate ID: EC50\_230619\_Plate1



**Friday 23-06-23**

**Split Cells**


1. Threw out BLLW (all dead)
2. Split JIRE 1/2
3. Split N2/RBL1 1/4

**Seed EC50**

- EC50 Plate ID: EC50\_230623\_Plate1

# 96 Well Plate Map

## EC50\_230623\_plate1




### Drug Dilutions

Drug	Dilution	Concentration	Units	Starting_Concentration
Cyclophosphamide	1	20	uM	20000
Cyclophosphamide	2	13.3	uM	13300
Cyclophosphamide	3	8.89	uM	8890
Cyclophosphamide	4	5.93	uM	5930
Cyclophosphamide	5	3.95	uM	3950
Cyclophosphamide	6	2.63	uM	2630
Cyclophosphamide	7	1.76	uM	1760
Cyclophosphamide	8	1.17	uM	1170
Cyclophosphamide	9	0.78	uM	780

Monday 26-06-23

### Collect EC50 Plate

- Plate ID: EC50\_230623\_Plate1



Break for writing MRes Report: 27/06/2023 - 24/07/2023

**JULY 2023**

**Monday 24-07-2023**

**Seeding Cells - RAMOS/FDC**

**Made New Media: AR-1**

1. Made new media: AR-1

Solution	ID code	Volume	% Total volume
Advanced RMPI	—	500 mL	78%
Glutamax	—	6.5 mL	20%
Pen-Strep	—	6.5 mL	1%
FBS	—	128 mL	1%

2. Seeded Feeders

- 2 vial of FDC cells thawed
- Resuspended in 12 mL of media
- 6 mL of FDC cell suspension were added to 2 separate T75 flasks (Seeding RBL1/BLLW on 25-07-2023)
- Remaining 12 mL of Media were added to T175 (Seeding N2 on 25-07-2023)

3. Seeded Ramos WILDseq Barcode Pool 2

- Thawed a vial of RAMOS pool 2 (previously barcoded DATE)
- Resuspended in 5 mL
- Spun down 90\* for 5min

- Resuspended in 12 mL of AR-1 media
- Seeded into T75

### Experimental Plans

1. Barcode BLLW/RBL1
2. Set up Drug pressure experiment for RAMOS barcoded cells
  - EC50's for cyclophosphamide and methotrexate need to be completed
  - Drugs to test (EC50):

Drug	Concentration	Replicates
DMSO	1%	3
Doxorubicin	200nM	3
Methotrexate	100nM	3
Vincristine	1nM	3

3. Do Side Population Experiment for N2 Barcoded Cells
  - Grow up sorted side pop to avoid issues with low cell numbers

## Tuesday 25-07-2023

### Seed Cells - N2-BC4/RBL1/BLLW/RAMOS-BC1

- Got vials from Chris
- Thawed, resuspended in 6 mL AR-1
- Spun down at 90\* for 5 min
- Discarded supernatant and resuspended pellet in 12 mL of AR-1
- Added RBL1/BLLW to T75's pre-seeded with FDC (24-07-2023)
- Added RAMOS-BC1 to new T75
  - RAMOS-BC2 from 24-07-2023 was dead
  - Was taken from -80, New RAMOS vial from LN
- Added N2 to T125 seeded with FDC (24-07-2023)

## Wednesday 26-07-2023

### Cell Culture

- Checked cells, all flasks looked healthy
- Added media
  - N2 - 20ml of AR-1
  - RAMOS - 10 mL of AR-1
  - RBL1 - 10ml of AR-1
  - BLLW - 10ml of AR-1

## Thursday 27-07-2023

### Split Cells - RAMOS

- Split Ramos 1:2
- Added 10ml of AR-1 to N2
- Plan to split RBL1/BLLW tomorrow

## Lucy Drug Randomizer

- Randomized drugs for Lucy Hare mouse experiment

Drug	Identification
Alectinib	A
Erda	B
Vehicle	C
Alectinib + Erda	D

## Friday 28-07-2023

Drug	Dilution	Concentration	Units	Starting Concentration
Cyclophosphamide	1	20	uM	20000
Cyclophosphamide	2	13.3	uM	13300
Cyclophosphamide	3	8.89	uM	8890
Cyclophosphamide	4	5.93	uM	5930
Cyclophosphamide	5	3.95	uM	3950
Cyclophosphamide	6	2.63	uM	2630
Cyclophosphamide	7	1.76	uM	1760
Cyclophosphamide	8	1.17	uM	1170
Cyclophosphamide	9	0.78	uM	780

## OCTOBER 2023

### Friday 06-10-2023

#### Cell Culture

##### Made New Media: AR-2

- Made new media: AR-2

Solution	ID code	Volume	% Total volume
Advanced RMPI	—	500 mL	78%
Glutamax	—	6.5 mL	20%
Pen-Strep	—	6.5 mL	1%
FBS	—	128 mL	1%

## Monday 09-10-2023

#### Cell Culture

##### Seeded A20

- Got A20 from Swetha box: 2nd rack in shared -80
- AR-2

### **Seeded FDC**

- Seeded feeders, taken from feeder box, irradiated by Chris
- Will seed N4 PDX tomorrow

## **Tuesday 10-10-2023**

### **Cell Culture**

#### **Seeded N4**

- Given by Chris -AR-2
- Seeded onto FDC cells from 09-10
- Future experiments:
  - Barcode PDX
  - EC50: Vincristine/Doxorubicin/Methotrexate/Cyclophosphamide (Active)

#### **Seeded RAMOS-BC**

- Seeded barcoded Ramos cells - Pool #4
  - From Chris
  - LN tank 1 - Rack 5 - box 2nd from Bottom
- AR-2
- Future Experiments:
  - Combination Drug pressure: Vin/Dox/Meth

#### **Checked A20**

- Look healthy
- ~40% confluency

## **Wednesday 11-10-2023**

### **Cell Culture**

#### **A20 Culture**

- ~70% confluent
- Split A20 in half
  - Half in 1xT75 and other half in 1xT175
  - T175 to be frozen down

#### **RAMOS-BC4 Culture**

- ~60% confluent
- Spun down at 400xg for 5 min and resuspended in 20ml AR-2

#### **N4 Culture**

- Look healthy, attached to feeders with some floaters
- Will add media tomorrow
- Plans:
  - Freeze down aliquots (2-3x10<sup>6</sup> per vial needed, ideally freeze down 3x vials)
  - Plan barcoding experiment

## Thursday 12-10-2023

### Cell Culture

#### Ramos-BC4

- Infected, will need to grab another vial
- Is the Stock infected? Will thaw 2 different pools this time

#### A20

- 1 x T75 + 1 x T175 both about ~60% confluent
- Media is slightly turbid, Jamie did a visual check and cells seem clean (A20 have a tendency towards slight turbidity - Jamie)
- Will add 10mL of AR-2 and observe till tomorrow
- If clean:
  - Freeze down 3 vials from T175
  - Set up EC50s with left over cells (Vin/Dox/Meth)
  - Split

#### N4

- Look healthy, media clean
- About ~40-50% confluent
- Will add 10mL of AR-2 and observe on Friday

## Friday 13-10-2023

### Cell Culture

#### A20

- Freeze down (10% DMSO in FBS)
- A20 count:  $2.91 \times 10^6$ /mL
  - 2mL + 3mL Media
- 3 vials frozen and put in LN Tank 6, Rack 5, Box 6

#### EC50\_231013

##### EC50 Plate ID: EC50\_231013\_Plate1

- Doxorubicin top/Vincristine bottom

Doxorubicin Stock: 10mM DMSO


Drug	Dilution	Concentration	Units	Starting_Concentration
Doxorubicin	1	300	nM	300000
Doxorubicin	2	200	nM	200000
Doxorubicin	3	133	nM	133000
Doxorubicin	4	88.9	nM	88900
Doxorubicin	5	59.3	nM	59300
Doxorubicin	6	39.5	nM	39500
Doxorubicin	7	26.3	nM	26300
Doxorubicin	8	17.6	nM	17600
Doxorubicin	9	11.7	nM	11700

Vincristine stock: 30 $\mu$ M

Drug	Dilution	Concentration	Units	Starting_Concentration
Vincristine	1	10	nM	10000
Vincristine	2	7.14	nM	7140
Vincristine	3	5.1	nM	5100
Vincristine	4	3.64	nM	3640
Vincristine	5	2.6	nM	2600
Vincristine	6	1.86	nM	1860
Vincristine	7	1.33	nM	1330
Vincristine	8	0.949	nM	949
Vincristine	9	0.678	nM	678

## 96 Well Plate Map

### EC50\_231013\_plate1



### EC50 Plate ID: EC50\_231013\_Plate2

- Methotrexate top/Cyclophosphamide bottom

Methotrexate Stock: 10mM DMSO

Drug	Dilution	Concentration	Units	Starting_Concentration
Methotrexate	1	500	nM	500000
Methotrexate	2	250	nM	250000
Methotrexate	3	125	nM	125000
Methotrexate	4	62.5	nM	62500
Methotrexate	5	31.2	nM	31200
Methotrexate	6	15.6	nM	15600
Methotrexate	7	7.81	nM	7810
Methotrexate	8	3.91	nM	3910
Methotrexate	9	1.95	nM	1950


Cyclophosphamide stock: 20 mM DMSO



Drug	Dilution	Concentration	Units	Starting_Concentration
Cyclophosphamide	1	20	uM	20000
Cyclophosphamide	2	14.3	uM	14300
Cyclophosphamide	3	10.2	uM	10200
Cyclophosphamide	4	7.29	uM	7290
Cyclophosphamide	5	5.21	uM	5210
Cyclophosphamide	6	3.72	uM	3720
Cyclophosphamide	7	2.66	uM	2660
Cyclophosphamide	8	1.9	uM	1900
Cyclophosphamide	9	1.36	uM	1360

## 96 Well Plate Map

### EC50\_231013\_plate2



Monday 16-10-2023

Cell Culture

A20

- Split 1/5
- AR-2

N4

- Split 1/5
- AR-2

Collect EC50\_231013

Plate ID: EC50\_231013\_Plate1

- Vincristine DNW




Figure 1: A20 Doxorubicin Results




Figure 2: A20 Methotrexate Results




Figure 3: A20 Cyclophosphamide Results

Plate ID: EC50\_231013\_Plate2

**Thursday 19-10-2023**

**Cell Culture**

**A20**

- Split 1/5
- AR-2

**N4**

- Split 1/5
- AR-2

**Friday 20-10-2023**

**Cell Culture**

**A20**

- Split 1/5
- AR-2

**N4**

- Split 1/5
- AR-2

## Sunday 22-10-2023

### Cell Culture

#### RAMOS Barcode

- LN Tank 6, Rack 5, Box 6
- Thawed Ramos-BC 1, 5, 6

## Monday 23-10-2023

### Cell Culture

#### N4

- ~40-50% confluent
- Added 10mL AR-2

#### A20

- ~60-70% confluent
- Added 10mL AR-2

#### Ramos-BC 1

- ~50% confluent
- Added 10mL AR-2

#### Ramos-BC 5

- ~50% confluent
- Added 10mL AR-2

#### Ramos-BC 6

- ~50% confluent
- Added 10mL AR-2

### Seed Feeders (FDC)

- Count:  $4.19 \times 10^5$
- Dilute to  $4.8 \times 10^4$
- Resuspended in 5mL AR-2, added 3.73 mL
- Seeded 1x 6-well plate
  - 2 mL/well
- Seeded 1x 96-well plate (15 wells)
  - 90  $\mu$ L/well

## Tuesday 23-10-2023

### Cell Culture

#### N4

- No split, maintained

## A20

- No split, maintained

## Ramos-BC 1

- ~60% confluent
- Split 1/3
- AR-2

## Ramos-BC 5

- ~60% confluent
- Split 1/3
- AR-2

## Ramos-BC 6

- ~60% confluent
- Split 1/3
- AR-2

## FDC

- Looked attached, healthy
- Media fine

## Transduce A20/N4 with WILDseq Library

1. Made 2mL cell suspension of  $4 \times 10^6$  cells/mL
  - A20:  $1.91 \times 10^6$  count
  - N4:  $1.3 \times 10^5$  count (took 2 mL of CS from flask)
2. Add 0.5mL of cell suspension to 4 wells of 12 well plate
3. Add virus mixture (listed below) to respective wells:

Well ID	Virus	Media	HEPES
A1	25 $\mu$ L	900 $\mu$ L	50 $\mu$ L
A2	10 $\mu$ L	930 $\mu$ L	50 $\mu$ L
A3	1 $\mu$ L	948 $\mu$ L	50 $\mu$ L
A4	0 $\mu$ L	950 $\mu$ L	50 $\mu$ L

4. Centrifuged cells at 600 x g for 1.5 hours at 32C
5. Resuspended cells and transferred to 6 well plate
  - N4 were put on feeder plate which was seeded 23-10-2023

## Wednesday 25-10-23

### Attempt 1: A20/N4 Barcoding - Media change

- Changed media
- AR-2

## Made New Media: AR-3

1. Made new media: AR-3

Solution	ID code	Volume	% Total volume
Advanced RMPI	—	500 mL	78%
Glutamax	—	6.5 mL	20%
Pen-Strep	—	6.5 mL	1%
FBS	—	128 mL	1%

## Thursday 26-10-23

### Attempt 1: A20/N4 Barcoding - Sort

- DNW
- No cells were positive
- Possible causes:
  - Virus volume too low
  - Spin speed too low
  - packaging envelope incorrect for mouse cells (GALV)

## Friday 27-10-23

### Cell Culture

#### A20

- Maintained
- No Split

#### N4

- Split 1/4
- AR-3

#### Ramos-BC 1

- ~60% confluent
- Split 1/4
- AR-3

#### Ramos-BC 5

- ~60% confluent
- Split 1/4
- AR-3

#### Ramos-BC 6

- ~60% confluent
- Split 1/4
- AR-3

## EC50\_231027

### Overview

- Treated A20 with Doxorubicin, Vincristine, Cyclophosphamide, Methotrexate
- Concentrations were the same as EC50\_231013
- Repeat for biological replicates
- Made an error for plate 1

### EC50 Protocol:

1. Create 12 ml cell suspension of  $2.5 \times 10^6$  cells
  - A20 Cell count:  $2.5 \times 10^6$
  - 8.33mL AR-3 added to 1.66mL of flask CS
2. Add 90  $\mu$ L of CS to respective wells of a 96-well plate
3. Dilute drug to starting concentration indicated in DMSO
4. Serial dilute into DMSO 8 times
5. Add 10uL of diluted drug to respective eppendorf containing 990uL media
6. Add 10 $\mu$ L of drugged media to respective well

### EC50 Plate ID: EC50\_231027\_Plate1

- Doxorubicin top/Vincristine bottom

Doxorubicin Stock: 10mM DMSO


Drug	Dilution	Concentration	Units	Starting_Concentration
Doxorubicin	1	300	nM	300000
Doxorubicin	2	200	nM	200000
Doxorubicin	3	133	nM	133000
Doxorubicin	4	88.9	nM	88900
Doxorubicin	5	59.3	nM	59300
Doxorubicin	6	39.5	nM	39500
Doxorubicin	7	26.3	nM	26300
Doxorubicin	8	17.6	nM	17600
Doxorubicin	9	11.7	nM	11700

Vincristine stock: 30 $\mu$ M DMSO

Drug	Dilution	Concentration	Units	Starting_Concentration
Vincristine	1	10	nM	10000
Vincristine	2	7.14	nM	7140
Vincristine	3	5.1	nM	5100
Vincristine	4	3.64	nM	3640
Vincristine	5	2.6	nM	2600
Vincristine	6	1.86	nM	1860
Vincristine	7	1.33	nM	1330
Vincristine	8	0.949	nM	949
Vincristine	9	0.678	nM	678

# 96 Well Plate Map

## EC50\_231013\_plate1



NOTE: Made an error by adding Methotrexate and Cyclophosphamide treatments to PLate 1 (as well as listed treatments) - Disregard results

### EC50 Plate ID: EC50\_231027\_Plate2

- Methotrexate top/Cyclophosphamide bottom

Methotrexate Stock: 10mM DMSO

Drug	Dilution	Concentration	Units	Starting_Concentration
Methotrexate	1	500	nM	500000
Methotrexate	2	250	nM	250000
Methotrexate	3	125	nM	125000
Methotrexate	4	62.5	nM	62500
Methotrexate	5	31.2	nM	31200
Methotrexate	6	15.6	nM	15600
Methotrexate	7	7.81	nM	7810
Methotrexate	8	3.91	nM	3910
Methotrexate	9	1.95	nM	1950


Cyclophosphamide stock: 20 mM DMSO



Drug	Dilution	Concentration	Units	Starting_Concentration
Cyclophosphamide	1	20	uM	20000
Cyclophosphamide	2	14.3	uM	14300
Cyclophosphamide	3	10.2	uM	10200
Cyclophosphamide	4	7.29	uM	7290
Cyclophosphamide	5	5.21	uM	5210
Cyclophosphamide	6	3.72	uM	3720
Cyclophosphamide	7	2.66	uM	2660
Cyclophosphamide	8	1.9	uM	1900
Cyclophosphamide	9	1.36	uM	1360

## 96 Well Plate Map

### EC50\_231013\_plate2



Monday 27-10-23

Cell Culture

A20

- Split 1/2
- AR-3

N4

- Split 1/2
- AR-3

Ramos-BC 1

- Split 1/2
- AR-3

## Ramos-BC 5

- Split 1/2
- AR-3

## Ramos-BC 6


- Split 1/2
- AR-3

## Seed Feeders

- Repeating barcoding of A20/N4
- Count:  $3.64 \times 10^5$
- Dilute to  $4.8 \times 10^4$
- Resuspended in 5mL AR-2, added 2.58 mL
- Seeded 1x 6-well plate
  - 2 mL/well
- Seeded 1x 96-well plate (15 wells)
  - 90  $\mu$ L/well

## Collect EC50\_231027

Plate ID: EC50\_231027\_Plate



Tuesday 31-10-2023

## Attempt 2: Transduce A20/N4 with WILDseq Library

1. Made 2mL cell suspension of  $4 \times 10^6$  cells/mL
  - A20:  $2.18 \times 10^6$  count
  - N4:  $3.06 \times 10^5$  count (took 2 mL of CS from flask)

2. Add 0.5mL of cell suspension to 4 wells of 12 well plate
3. Add virus mixture (listed below) to respective wells:

Well ID	Virus	Media	HEPES
A1	250 $\mu$ L	225 $\mu$ L	25 $\mu$ L
A2	150 $\mu$ L	150 $\mu$ L	25 $\mu$ L
A3	50 $\mu$ L	50 $\mu$ L	25 $\mu$ L
A4	0 $\mu$ L	0 $\mu$ L	25 $\mu$ L

4. Centrifuged cells at 1000 x g for 1.5 hours at 0C

- Mistakenly set centrifuge to 0C

5. Resuspended cells and transferred to 6 well plate

- N4 were put on feeder plate which was seeded 23-10-2023

## Drug Pressure Experiment 1: RAMOS Day 1

### Overview

- Designing combination treatment to reflect common BL treatments in clinic
- Plan to analyse differences in barcode representation between baseline, DMSO, and treatment-resistant
- Drug concentrations:

Drug	Treatment Day	Concentration	Stock
Doxorubicin	Day 2	80nM	100 $\mu$ M working stock
Vincristine	Day 1	10nM	30 $\mu$ M
Methotrexate	Day 1	20nM	100 $\mu$ M working stock
Cyclophosphamide	Day 2/3/4	20 $\mu$ M	20 $\mu$ M working stock


- Plan:
  1. Collect sample at baseline
  2. Complete 4 day treatment course of RAMOS-BC cells (with corresponding DMSO controls)
  3. Allow treatment(T) cells to recover and collect sample
  4. Complete 4 day treatment course again and collect sample
- Possible future experiments:
  - EC50's of T vs WT cells (are T RAMOS cells drug resistant?)

## Day 1

### Seeded RAMOS BC

- Seeded 3x 6 well plates with RAMOS-BC 1, 5, 6 respectively

## 96 Well Plate Map



- Seeded  $0.25 \times 10^6$  cells/well in 3mL AR-3
- Treatment wells

Drug	Treatment Day	Concentration	Stock	Stock Volume	Media Volume
DMSO	Day 1	-	-	14.4 $\mu$ L	27 mL AR-3
Vincristine	Day 1	10nM	30 $\mu$ M	9 $\mu$ L	27 mL AR-3
Methotrexate	Day 1	20nM	100 $\mu$ M working stock	5.4 $\mu$ L of WS	27 mL AR-3

### Collected RAMOS Baseline

- Took 1mL CS from respective flask and transferred to eppendorf
- Spun down 2500 x g at 4C
- Removed Media
- Snap froze in LN
- Put in -80 CM Box 1

**NOVEMBER 2023**

**Wednesday 01-11-2023**

### Drug Pressure Experiment 1: RAMOS Day 2

- Changed media
- Spun down at 600 x g for 5 min
- Resuspended in Day 2 Treatment:

Drug	Treatment Day	Concentration	Stock	Stock Volume	Media Volume
DMSO	Day 1	-	-	48.60 $\mu$ L	27 mL AR-3
Doxorubicin	Day 1	80nM	100 $\mu$ M working stock	21.60 $\mu$ L of working stock	27 mL AR-3
Cyclophosphamide	Day 1	20 $\mu$ M	20mM	27 $\mu$ L	27 mL AR-3

## Attempt 2: Barcoding A20/N4

- Sanity check
- Ran samples on E6 flow
- A20 negative
- Suggests envelope protein is not compatible with mouse cells
- N4 dead
- Re-make WILDseq virus with VSVG packaging envelope

## Thursday 02-11-2023

### Cell Culture

#### A20

- Split 1/2
- AR-3

#### N4

- Split 1/2
- AR-3

#### Ramos-BC 1

- Split 1/4
- AR-3

#### Ramos-BC 5

- Split 1/4
- AR-3

#### Ramos-BC 6

- Split 1/4
- AR-3

### Seeded Feeders

- 1 vial
- Split into 2 T75 flasks

## Drug Pressure Experiment 1: RAMOS Day 3

- Changed media
- Spun down at 600 x g for 5 min
- Resuspended in Day 3 Treatment:

Drug	Treatment Day	Concentration	Stock	Stock Volume	Media Volume
DMSO	Day 1	-	-	27 $\mu$ L	27 mL AR-3
Cyclophosphamide	Day 1	20 $\mu$ M	20mM	27 $\mu$ L	27 mL AR-3

## Friday 03-11-2023

### Cell Culture

#### A20

- Split 1/2
- AR-3

#### N4

- Split into two new feeder flasks
  - Seeded 02-11-2023
  - N4.1 / N4.2
  - 1 flask to transduce, 1 flask to freeze down
- AR-3

#### Ramos-BC 1

- Add 10mL AR-3

#### Ramos-BC 5

- Add 10mL AR-3

#### Ramos-BC 6

- Add 10mL AR-3

### Drug Pressure Experiment 1: RAMOS Day 4

- Treatment cells too sick to continue
- Changed media
- Spun down at 600 x g for 5 min
- Resuspended both DMSO and Treatment cells in fresh AR-3
- Allow to grow

## Monday 06-11-2023

### Cell Culture

#### A20

- Split 1/2
- AR-3

#### N4.1

- Healthy/growing
- Leave till tomorrow

## N4.2

- Healthy/growing
- Leave till tomorrow

## Ramos-BC 1

- Froze down
- 3x vials
- FBS +10% DMSO

## Ramos-BC 5

- Froze down
- 3x vials
- FBS +10% DMSO

## Ramos-BC 6

- Froze down
- 3x vials
- FBS +10% DMSO

## Seed Feeders

- For WILDseq transduction of N4 (tuesday?)
- 1x 4 wells of 6-well plate, 15wells of 96-well plate
- Dilute in 8ml
- 2mL/well of 6well
- 90 $\mu$ L/well of 96well plate

## Drug Pressure Experiment 1: RAMOS Day 7

- Split DMSO cells 1/6
- Add 1mL AR-3 to treatment cells
  - treatment cells still sick, allow to grow for a few more days
  - Perhaps transfer to smaller well?

**Tuesday 07/11/2023**

## Made New Media: AR-4

1. Made new media: AR-4

Solution	ID code	Volume	% Total volume
Advanced RMPI	—	500 mL	78%
Glutamax	—	6.5 mL	20%
Pen-Strep	—	6.5 mL	1%
FBS	—	128 mL	1%

## Cell Culture

### A20

- Split 1/6
- AR-4

#### N4.1

- Healthy/growing
- Leave till tomorrow

#### N4.2

- Froze down
- 2x vials
- FBS +10% DMSO

### Drug Pressure Experiment 1: RAMOS Day 5

- Split D 1/6

#### Attempt 3: Transduce N4 with WILDseq Library

1. Made 2mL cell suspension of  $4 \times 10^6$  cells/mL
  - N4: took 10ml of 20ml CS, spun down and resuspended in 2mL AR-4
2. Add 0.5mL of cell suspension to 4 wells of 12 well plate
3. Add virus mixture (listed below) to respective wells:

Well ID	Virus	Media	HEPES
A1	250 $\mu$ L	225 $\mu$ L	25 $\mu$ L
A2	150 $\mu$ L	150 $\mu$ L	25 $\mu$ L
A3	50 $\mu$ L	50 $\mu$ L	25 $\mu$ L
A4	0 $\mu$ L	0 $\mu$ L	25 $\mu$ L

4. Centrifuged cells at 1000 x g for 1.5 hours at 32C
5. Resuspended cells and transferred to 6 well plate
  - N4 were put on feeder plate which was seeded 06-11-2023

### Wednesday 08/11/2023

#### Attempt 3: Barcoding N4

- Changed Media
- Cells look healthy ## Drug Pressure Experiment 1: RAMOS Day 6
- Checked treatment cells with trypan blue
- Everything was dead
- Ended experiment: will reconsider approach
  - Treat Ramos with individual agents as well as combination
  - Treat for 3 days and recover or 12 days and recover(?)

### Thursday 09/11/2023

#### Attempt 3: Barcoding N4

- Booked sort for tomorrow
- Changed media
- Sanity check on E6 Fortessa (very low MOI) ## Cell Culture



## A20

- Split 5mL CS into 15mL of AR-4 ### N4
- Split 5mL CS into 30mL of AR-4

## Friday 10/11/2023

### Attempt 3: Barcoding N4

- Sorted N4
- 1000 cells into wells A1-A3 of 96-well plate seeded with feeders
  - Seeded 06/11/2023
- 500 cells into wells B1-B3
  - Seeded 06/11/2023

## Cell Culture

### RAMOS-BC1

- Thawed vial of RAMOS-BC1 into 15mL of AR-4 ### N4
- Threw away ### A20
- Threw away

## Monday 13/11/2023



### Start RAMOS-Drug Pressure v2


1. Count RAMOS-BC1 CS
  - CS:  $1 \times 10^6$  cells/mL
2. Add 122.8 $\mu$ L CS per well
3. Make drug dilutions according to following table:
  - Make working stock of Doxorubicin and Methotrexate
  - Dilute 10mM stock to 100 $\mu$ M
    - 10 $\mu$ L stock in 990  $\mu$ L DMSO

DMSO: 16.67  $\mu$ L DMSO in 10mL Combo: Added doxorubicin, vincristine, methotrexate, and cyclophosphamide to 10mL at same volumes indicated below

Drug	Concentration	Stock	Dilution
Doxorubicin	80nM	100 $\mu$ M working stock (10mM Stock)	5 $\mu$ L WS in 10mL
Vincristine	10nM	30 $\mu$ M	1.67 $\mu$ L stock in 10mL
Methotrexate	20nM	100 $\mu$ M working stock (10mM Stock)	5 $\mu$ L WS in 10mL
Cyclophosphamide	20 $\mu$ M	20 $\mu$ M working stock	5 $\mu$ L stock in 10mL

4. Add 3 mL of drugged media to respective wells





## Cell Culture

### N4-BC Attempt 3

- Add media to wells with cells
- Lots of death, some clusters of live cells ### RAMOS-BC1
- Split 4mL CS into 20mL AR-4

### PIL COURSE

**Thursday 13/11/2023**

### RAMOS-DPv2

Change media: 1. Transfer into 15ml Falcon 2. Spun down @ 500xg for 6min 3. Resuspended in 1mL 4. Transferred into new 12-well plate

**Friday 14/11/2023**

## Cell Culture

**INFECTION** - RAMOS-BC1 and RAMOS-DPv2 infected - Threw out cultures and AR-4 - N4 seemingly unaffected